

All-in-One HD Video Magnifier with Selective Text-to-Speech

DaVinci is a lightweight desktop video magnifier (CCTV) that can be easily transported from room to room. Experience Full HD, selective text-to-speech (OCR) and a high-performance 3-in-1 camera. High definition color and contrast provides beautiful crystal-clear images and vibrant colors. Experience the joy of reading with our new text-to-speech feature. DaVinci will read any printed text aloud with the push of a button.

DaVinci Features:

- Full HD 1080p camera
- 24" high resolution HD LCD for best image quality
- Auto focus camera for self-viewing, reading or distance viewing (3-in-1 camera)
- Selective Text-to-Speech (OCR)
- Magnification up to 77x
- Easy-to-use console with line markers
- Reading zone preview with motion detection
- 28 viewing modes to optimize contrast and brightness
- Slide mechanism provides flexibility for various camera arm positions
- LCD screen easily pivots vertically
- Memory setting for each camera position
- Download saved images to a flash drive
- DaVinci's monitor can also be utilized as a display for your computer or iPad (additional hardware required)
- Set up is easy, just plug in and begin
- 2 year warranty



“The magnification and text-to-speech features make reading enjoyable again!”

– Jerry Winkler



Camera rotates 360°

See near, far and everything in between.

Easily pivot DaVinci's screen to provide a comfortable viewing position.

enhanced
vision®
A VISPERO™ BRAND

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

All-in-One HD Video Magnifier with Selective Text-to-Speech



DaVinci connected to iPad



- Full HD 1080p camera displays crystal clear images in vibrant color and contrast
- High resolution LCD produces maximum picture detail for a clear, bold display
- Large field of view – see more on the screen
- Let DaVinci read your favorite article or book aloud
- Place printed text under the camera and press capture button
- DaVinci begins reading what's on the screen in a few seconds
- Many different languages are available
- See Near, Far and Mirror image with DaVinci's 3-in-1 camera
- Apply make-up, shave, read, write, view presentations and whiteboards
- Perfect for work, home or school

Local Dealer:

To schedule a no obligation demonstration, please call:

enhanced
vision®
A VISPERO™ BRAND

(888) 811-3161

(714) 465-3400

www.enhancedvision.com