

Portable HD Desktop Magnifier On The Go!

Amigo HD displays amazing crystal clear HD images with vibrant colors, at any distance, on a large 7 inch widescreen LCD. So lightweight and simple to use, you will never want to be without it. Place it directly on your reading material or hold it up to view an item at a distance. Adjust magnification and contrast with an easy press of a button. Read price tags, menus, labels, street signs and more, anywhere!

Amigo HD Features:

- HD camera produces vibrant colors and crystal clear images
- 7" High resolution LCD with large field of view
- Auto focus & one push manual focus for viewing objects at any distance
- Adjustable magnification 1.4x to 25x
- Easy-to-use large tactile buttons
- Adjustable brightness with lights on/off
- 28 available color modes
- Adjustable folding base
- Built-in rechargeable lithium ion battery with 3 to 5 hours battery life (depending on the mode and brightness level)
- Low battery display indicator
- Image capture & save with capability to magnify, pan, change viewing modes & upload to PC
- Video output to a monitor or television
- Clock feature
- Weighs only 1.3 pounds
- Buttons available in Color or Black & White
- Carrying case included
- 2 year warranty



“My new Amigo HD is the perfect companion when I am on the go!”

– Barbara Johnson



**enhanced
vision**
(888) 811-3161
(714) 465-3400
www.enhancedvision.com

Portable HD Desktop Magnifier On The Go!

Amigo HD is the perfect choice for reading mail, magazines, newspapers, prescriptions, books, viewing family photos, seeing objects at varying distances and saving captured images for later viewing.



Amigo HD's lightweight portable design makes it the perfect companion when you are at home or on the go.



Amigo HD offers two desktop hands free positions, allowing you to easily read, write, pay bills, file nails and work on hobbies.



Local Dealer:

To schedule a no obligation demonstration, please call:

enhanced
vision

(888) 811-3161

(714) 465-3400

www.enhancedvision.com