

transformer **HD**

with Built-in Wi-Fi

PORTABLE ELECTRONIC
MAGNIFIER WITH BUILT-IN WI-FI
AND OPTIONAL OCR



A small but powerful electronic magnifier with self-view, up-close, and distance viewing capabilities, the Transformer HD offers built-in Wi-Fi, as well as text-to-speech on the OCR model. Easily transport this lightweight device from home to the office, school, or anywhere you want to go.

More Choices. More Solutions. You'll see.

enhanced
vision[®]
A VISPERO[®] BRAND

transformer **HD**

with **Built-in Wi-Fi**

High-performance magnification on the go



CCTV Model Features:

- Full HD 1080p auto-focus 3-in-1 camera for self-viewing, reading, or distance viewing
- Connects to your monitor, computer, or tablet via HDMI, USB 3.0, or Wi-Fi
- 4 desktop positions available
- Magnification up to 30x on 24" LCD
- 32 viewing modes to optimize contrast and brightness
- Work in full- or split-screen mode
- Removable legs for suction base option
- Record and save videos
- Built-in LED lighting
- 3-hour minimum removable battery with internal charger (spare battery and external charger included)
- Memory setting for each camera position

CCTV with OCR Model Features:

- Optional 13-megapixel camera on the OCR model for accurate, full-page text-to-speech
- Adjustable magnification and color viewing modes in OCR
- Reading preview lets you read the entire page, or use your cursor to select where to begin
- Multiple language recognition with male and female premium voices
- Save and recall documents, pictures, and export files to your computer

enhanced
vision
A VISPERO® BRAND

Local Dealer:



Take Transformer HD anywhere with its convenient folding design. Carrying case included.

MKT-0126-00, Rev. 08/2023

2-year warranty included

888-811-3161 | 714-465-3400 | www.enhancedvision.com