



Compact • Portable • Lightweight

## The only Hand-Held Magnifier you will ever need!

Pebble's small lightweight design makes it the perfect companion whether you're travelling or at home. Easily read labels, prescriptions, price tags, menus, bus schedules and so much more with this ultra-portable video magnifier.

Carry Pebble in your pocket, purse or strap it onto your belt, the only time you'll know it's there is when you need it!

### Pebble Features:

- Adjustable magnification from 2x to 10x
- Lightweight design
- 3.5" or 4.3" high resolution LCD
- Large viewing area in a small package
- Adjustable brightness
- Easy-to-use tactile controls
- 28 available viewing modes
- Freeze image feature with capability to magnify and change mode
- Foldable handle with comfort grip
- Easy writing capability
- Auto focus
- Over 2 hours battery life (rechargeable/replaceable batteries included)
- On/Off light option for reduced glare
- 2 year warranty



“ I really enjoy the larger screen on my Pebble 4.3, it allows me to see so much more and it's nice to carry one magnifier! ”

– Vicki Williams,  
Bank Teller



**enhanced**  
**vision**®

**(888) 811-3161**

**(714) 465-3400**

[www.enhancedvision.com](http://www.enhancedvision.com)



or



Available in 3.5" or larger 4.3" LCD

	Pebble	Standard Optical Magnifiers
• Adjustable Magnification	✓	—
• 28 Available Viewing Modes	✓	—
• Freeze Frame Feature	✓	—
• Foldable Handle	✓	—
• Fully Portable	✓	✓
• Large Field of View	✓	—
• Writing Position	✓	—
• Adjustable Brightness	✓	—



To schedule a no obligation demonstration, please call  
(888) 811-3161

**enhanced**  
**vision®**

**(888) 811-3161**

**(714) 465-3400**

[www.enhancedvision.com](http://www.enhancedvision.com)

**Local Dealer:**