

A Full Page Smart Solution for Easy Text Reading!

Finally, a Full Page reader that makes reading enjoyable again! Easy-to-use, Smart Reader boasts a lightweight portable sleek design with a user friendly console. Position your document and press a button, in seconds Smart Reader begins reading aloud. Connect Smart Reader to a TV or monitor and enlarge your captured image to follow along as it reads. For a true portable experience, attach Smart Reader to our new **13.3**" **Portable Battery Operated Monitor** and enjoy up to **3.5 hours** of reading anywhere! Perfect for work, home, school, vacation and on the go.



Smart Reader Features:

- Dual Function High Resolution System; Full Page reader only or connect to 13.3" Portable Monitor, TV or desktop monitor & read along
- User friendly Control Console with tactile buttons
- Capture an image, magnify and read along comfortably in page, line or word format
- 28 color modes available to optimize contrast and brightness
- Book mode for saving multiple pages
- Table mode allows you to easily read documents with tables
- Audio and visual prompts

- Volume and Speed controls
- Capture, save and recall modes
- Natural sounding male or female voice
- Over 20 languages available
- Clock with voice prompt on demand
- Attach headphone for privacy
- Weighs under 5 lbs.
- 2 year warranty

enhanced vision[®] (888) 811-3161 (714) 465-3400 www.enhancedvision.com



Portable • Easy-to-Use • Sleek Design

A Full Page Smart Solution for Easy Text Reading!

Additional Features & Benefits:

- Smart Reader is a Full Page OCR system that is ideal for reading magazines, books, newspapers, recipes and any other printed material. Simply position your document and press a button, in a few seconds Smart Reader begins to read aloud.
- Connect Smart Reader to 13.3" Portable Monitor, or any desktop monitor or TV (HDMI & DVI) and view the enlarged image in capture mode. Read along and listen as Smart Reader reads aloud.
- A convenient built-in handle and lightweight design allows for easy transport anywhere you go.



Designed and Assembled in the USA



Microphone Accessory (Sold Separately)

Local Dealer:

Smart Reader is a stand alone OCR system that can read back and display an enlarged image in capture mode by connecting to a monitor or a large TV.

Smert Reader is a stand storm CCR system mode by connecting In the capture mode by connecting In the capture mode by connecting

Monitor Sold Separately

To schedule a no obligation demonstration, please call (888) 811-3161 (714) 465-3400 www.enhancedvision.com