

DaVinci connects to iPad & iPhone

DaVinci™ 

The Revolutionary All-in-One HD CCTV

with Text-to-Speech(OCR) and 3 camera viewing positions

See near, far and everything in between.

- **Full HD 1080p camera** - with full HD you will experience amazing picture detail
- **Auto focus camera** - for reading, distance viewing or self-viewing
- **Reads aloud with Nuance® software** - read any printed text aloud with the push of a button
- **24" high resolution HD LCD** - for crystal clear pictures and vibrant colors
- **28 adjustable viewing modes** allow you to personalize DaVinci - perfect for reading magazines, newspapers, etc.
- **Magnification** - up to 77x



Here is all you need:

Apple iPad Dock Connector to VGA (female) Part # MC552ZM/B or if you are using 4th generation iPad you need a cable called Lighting to VGA (female) MD825ZM/A. Also required is a VGA cable (male to male). One end connects to the Apple cable and the other connects to the VGA connector at the back of the DaVinci or Merlin elite monitor.



CALL NOW FOR MORE INFORMATION AND YOUR FREE DEMONSTRATION

888-811-3161

www.enhancedvision.com