

3-IN-1 CCTV FULL-HD DESKTOP VIDEO MAGNIFIER WITH OCR





DaVinci Pro is a high-performance desktop video magnifier (CCTV) featuring a full HD 1080p 3-in-1 camera and full-page text-to-speech (OCR). A lightweight design allows you to easily transport it from room to room.






DaVinci Pro HD OCR's 3-in-1 full HD camera makes it the perfect electronic magnifier for any task. See near, far, or even view yourself—apply make-up, shave, read, write, view presentations and whiteboards, and so much more! With the simple press of a button, toggle between CCTV and OCR (text-to-speech) modes for ultimate ease of use.

CCTV Mode Features:

- Auto-focus 3-in-1 1080p camera for close viewing, distance viewing, or self-view
- 24" high-resolution HD LCD screen with large field of view, pivots for best position
- Magnification up to 77x
- Slide mechanism provides flexibility for various camera arm positions
- Single-line text reading format to reduce eye fatigue
- 28 viewing modes to optimize contrast and brightness
- Computer compatible—toggle between CCTV and other devices
- Includes popular games like Sodoku, Soiltaire, and Tic-Tac-Toe



MKT-0123-00, Rev. 08/2023

OCR (Text-to-Speech) Mode Features:

- 13-megapixel camera for accurate full-page OCR—let DaVinci Pro HD read your favorite book or article aloud
- Reading preview lets you read the entire page or navigate through paragraphs and view selected text
- Easy-to-use console with simple and advanced modes
- Male and female premium voices
- Easily switch between languages
- Save documents, pictures or books and export files to your PC
- Setup is simple—just plug in and begin

enhanced Vision

Local D	ealer:
---------	--------

2-year warranty included