

transformer

VGA/USB Compatible • Portable • Lightweight

The next generation of portable electronic magnifiers:

Transformer is the most flexible and portable solution for reading, writing and viewing magnified images at any distance. Compatible with popular magnification software and the latest Windows* operating systems as well as Mac OS X*, Transformer's unique compact design is perfect for school, work, or home.

Transformer Features:

- Weighs less than three pounds
- 2.4x to 30x magnification (based on 17" screen)
- Camera rotates 330° for reading, distance & self-viewing modes
- Auto install software/plug & play
- Compatible with popular magnification software & the latest Windows* operating systems as well as Mac OS X*
- USB connectivity for laptop/computer
- VGA connectivity direct to monitor
- Capture and store images (when connected to computer)
- 28 custom color select modes
- Built in LED lighting
- Battery operated – up to 4 hours
- Includes soft carrying case & cloth sleeve
- Memory setting for each camera position
- 2 year warranty



Transformer connects to monitor, laptop or desktop computer!



enhanced
vision®

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

*Windows is a registered trademark of Microsoft Corporation. All other brand names are trademarks of their respective owners.

*Mac OS X is a registered trademark of Apple Inc. All other brand names are trademarks of their respective owners.

transformer

VGA/USB Compatible • Portable • Lightweight



Transformer's sleek compact design and traveling case is ideal for easy transport.



To schedule a no obligation demonstration, please call (888) 811-3161.

Local Dealer:



Connect to a PC laptop and Mac OS X*



Connect to a desktop computer or a monitor!



Transformer VGA also available with a battery operated portable 13.3" tilting monitor.

enhanced
vision®

(888) 811-3161

(714) 465-3400

www.enhancedvision.com

08/14

Designed and Assembled in the USA