

### **PORTABLE HIGH-DEFINITION** VIDEO MAGNIFIER ON THE GO

Featuring Dynamic Contrast<sup>™</sup> for sharper, clearer more defined images and text

## **NEW ENHANCED DESIGN AND FEATURES!**



The Amigo is a lightweight, portable video magnifier with a generous 8-inch screen and colorful tactile buttonseasily magnify text and objects up close, or even at a short distance.

ous shape by January 196 as simple: change it but do sar design exec Don Kopk pciety of Automotive Engi

n October 1966

More Choices, More Solutions, You'll see.



# amigo

Lightweight and easy to use, simply place Amigo directly on your reading material, or hold it up to view objects at a short distance. Adjust magnification and contrast with the simple press of a button. Read price tags, menus, labels, street signs, and more—anywhere!

### **Amigo Features:**

- 8" high-resolution LCD screen
- 2 high-definition cameras for up-close and short distance viewing
- Live image panning—pan through a live image without moving your material
- Built-in rechargeable battery provides up to 3.5 hours of continuous use
- Built-in reading stand positioned at a comfortable viewing angle
- Adjustable magnification up to 30x
- Line marker and masking features
- 24 available color modes
- Image capture with save capability
- Video output over USB-C port
- Easy-to-use, large tactile buttons
- Carrying case included



MKT-0129-00, Rev. 08/2023

#### All-New Visual Display Settings:

- Dynamic Contrast: provides a sharper, clearer, and more defined image, making undetectable text visible again
- Edge Effects: makes the edges of objects in an image stand out
- Color Blindness Filter: assists with color differentiation in real time
- Light Sensitivity Filter: mutes the lighting of very bright colors



Local Dealer:

2-year warranty included