

merlin **FULL HD**  
ultra —

FULL-HD DESKTOP VIDEO MAGNIFIER



*The affordable Merlin ultra offers full high-definition color and contrast, resulting in sharp, crystal-clear images and vibrant colors. Simply the greatest value and best picture quality available in HD desktop magnification—there's even more to see with Merlin ultra!*


More Choices. More Solutions. You'll see.

**enhanced**  
**vision**®  
A VISPERO® BRAND

The affordable Merlin ultra offers full high-definition color and contrast, resulting in sharp, crystal-clear images and vibrant color accuracy. Merlin ultra's Full HD camera allows for a wide field of view, displaying more text on the screen in amazing detail. Easily read books and magazines, read and write letters, manage your financial records, view color photos, do crossword puzzles, and even pursue a favorite hobby.

### Merlin HD ultra Features:

- Full HD camera: 1920 x 1080 resolution in 22" & 24" LCD
- High-resolution HD LCD screen with a wide field of view for best image quality
- Adjustable magnification up to 73x (depending on screen size)
- Built-in XY table with a user-friendly lock mechanism for smooth, continuous reading
- Screen easily pivots horizontally and vertically to provide the most comfortable viewing position
- Simple, easy-to-use tactile controls
- Adjustable color select modes for increased contrast personalization
- Horizontal line markers and horizontal windowing to reduce eye fatigue
- LED lighting for truer picture quality
- Computer compatible—toggle between CCTV and computer (additional hardware may be required)
- Setup is easy—just plug in and begin



*Merlin HD ultra's XY table allows you to move your material horizontally and vertically for smooth, continuous viewing.*

**enhanced**  
**vision**  
A VISPERO® BRAND

**Local Dealer:**