

## ALL-IN-ONE CCTV FULL HD DESKTOP VIDEO MAGNIFIER WITH OCR





Merlin elite Pro is the perfect solution for anyone looking for a powerful, all-in-one electronic magnifier. View crystal-clear images in vibrant colors on the highresolution LCD, or listen as Merlin elite Pro reads your favorite book or article aloud.





Merlin elite Pro is a high-performance desktop video magnifier (CCTV), featuring a full HD 1080p camera and full-page text-to-speech (OCR). Its high-definition camera displays crystal-clear images in vibrant color and contrast, while the high-resolution LCD produces maximum detail on a clear, bold display. Easily toggle between live image CCTV mode and full-page OCR mode with the push of a button—listen along as Merlin elite Pro reads your favorite magazine or book aloud!

## **CCTV Mode Features:**

- Full HD 1080p auto-focus camera
- 24" high resolution HD LCD for superior image quality and a large field of view see more on the screen
- Magnification up to 70x
- 32 viewing modes to optimize contrast and brightness
- Screen easily pivots horizontally and vertically to provide the most comfortable viewing position
- Low-profile XY table with a user-friendly lock mechanism allows you to easily move your material horizontally and vertically
- Computer compatible—toggle between Merlin elite Pro's CCTV and your computer, connects to tablets & iPad



MKT-0125-00, Rev. 08/2023

## **OCR Mode Features:**

- 13-megapixel camera for accurate full-page text-to-speech (OCR)
- 3 OCR reading formats (full-page text, fullpage picture, and single-line text) for an optimum reading experience
- Easy-to-use console with simple and advanced modes
- Multiple language recognition
- Male and female premium voices
- Save documents, pictures or books and export files to your PC
- Setup is easy—just plug in and begin



3-year warranty included