

merlin **FULL HD** -mini

3-IN-1 PORTABLE, FOLDABLE  
DESKTOP VIDEO MAGNIFIER



*Merlin mini offers the same reading comfort, ease of use, and functionality as a desktop video magnifier but in a compact and foldable design. A 3-in-1 rotating camera offers up close, distance, and self-viewing magnification modes. Connect to external devices with USB-C connection.*

More Choices. More Solutions. You'll see.

**enhanced**  
**vision**<sup>®</sup>  
A VISPERO<sup>®</sup> BRAND

*Crisp, clear magnification wherever you are*

Merlin mini is the ideal solution for people who need magnification in different locations. Read text and view objects at home, in the office, classroom, or on the go. Its foldable design makes it extremely easy to take Merlin mini with you, and store it when not in use.

## Merlin mini Features:

- 3-in-1 rotating camera allows you to read text and view objects up close, at a distance, and even view yourself
- Crisp, clear full HD image—1080p with auto-focus
- Generous 15.6" or 17.3" screen that easily folds away
- Portable—easy to transport and use in different locations
- Unique zoom and pan features
- 3 selectable magnification ranges to suit your vision needs
- Read full pages without repositioning your material thanks to a moveable camera
- Connect iPhone® and Android™ devices for an enlarged screen image\*
- Connect to your TV screen via HDMI for an even larger image
- Up to 5-hour battery life



*Lightweight, easy to carry*



*Easy to fold and store*



*Connect to iPhone® and Android™*

*(\*refer to user manual for device compatibility)*

**Local Dealer:**